

# DIG Easy Hyb

Hybridization solution for nucleic acid blots with digoxigenin-labeled probes

Cat. No. 1 603 558 500 ml

Version 2, March 2003

Store at 15 - 25° C

## 1. Product overview

### Contents

Bottle	Label	Content
1	DIG Easy Hyb	500 ml, ready to use solution DNase and RNase free

### Product description

DIG Easy Hyb is non-toxic and does not contain formamide, yet the hybridization temperature should be calculated with the same equation that is used for buffers containing formamide, 50%.

### Application

DIG Easy Hyb can be applied for all types of nucleic acid blot hybridizations. It is especially designed for the use with digoxigenin (DIG)-labeled nucleic acid probes to targets bound to nylon membrane.

### Prehybridization

Prehybridization with DIG Easy Hyb is performed for 15-30 min at the appropriate hybridization temperature.

### Hybridization

DIG Easy Hyb can drastically reduce hybridization times to only 1-6 hours, depending on the type of hybridization. Only for high sensitivity requirements, overnight (12-16 h) hybridization is recommended.

Application	Recommended hybridization time
DNA fingerprinting (multiple locus probes)	2-4 hours
DNA fingerprinting (single locus probes)	overnight
Colony/plaque hybridization	>4 hours
Single-copy gene detection in human genomic blots	overnight
RNA:RNA hybridization	6 h - overnight
Oligonucleotide probes	1-6 hours

### Storage/ stability

The unopened bottle is stable at 15-25°C through the control date printed on the label.

### Storage of DIG-labeled probes in DIG Easy Hyb

Probes labeled with digoxigenin can be stored in DIG Easy Hyb at -15 to -25°C and be re-used several times, after denaturing at 68° C prior to use, **do not boil**.  
**Note:** Storage and reuse of RNA probes in DIG Easy Hyb is generally not recommended and depends on probe inherent stability and RNase-free handling conditions.

## 2. Procedures and required materials

### 2.1 Before you begin

#### Additional equipment required

Hybridization can be performed in temperature resistant, sealable

- Hybridization bags\* (Cat. No. 1 666 649)
- plastic or glass boxes
- petri dishes
- roller bottles

### 2.2 DNA : DNA hybridization

#### Hybridization temperature

The appropriate hybridization temperature is calculated according to GC content according to the following equation:

$$T_m = 49.82 + 0.41 (\% G + C) - (600/l)$$

[l = length of hybrid in base pairs]

$$T_{opt.} = T_m - (20 \text{ to } 25^\circ\text{C})$$

(The given numbers of the equation were calculated according to a standard equation for hybridization solutions containing formamide, 50%.)

The actual hybridization temperature  $T_{opt.}$  for hybridization with DIG Easy Hyb is 20-25° C below the calculated  $T_m$  value.  $T_{opt.}$  can be regarded as a stringent hybridization temperature allows up to 18 % mismatches between probe and target. When the degree of homology of your probe to template is less than 80%, you should lower  $T_{opt.}$  accordingly (approx. 1.4°C below  $T_m$  per 1 % mismatch) and also adjust the stringent washing steps accordingly (*i.e.* increase SSC concentration and lower washing temperature).

#### Example

For hybridization of human genomic DNA with a 100% homologous probe use 37-42°C, depending on the GC contents of the probe.

#### Procedure

In the following table the procedure for a DNA:DNA hybridization is described.

**Note:** Do not use open trays when working with DIG Easy Hyb buffer.

Step	Action
1	Pre-heat appropriate volume of DIG Easy Hyb (approx. 20 ml/100 cm <sup>2</sup> ) to hybridization temperature.
2	Incubate the blot for 15-30 min with gentle agitation. <b>Note:</b> The membrane should be well immersed and covered with DIG Easy Hyb.
3	Denature DIG-labeled DNA probe (5-25 ng/ml hybridization solution) by boiling for 5 min and rapidly cooling on ice-water.
4	Add to pre-heated DIG Easy Hyb (at least 3.5 ml/100 cm <sup>2</sup> membrane) and mix well but avoid foaming (bubbles may lead to background).
5	Pour off prehybridization solution and immediately add probe/DIG Easy Hyb mixture to membrane. <b>Note:</b> Do not add concentrated probe directly to avoid localized background.
6	Incubate with gentle agitation for at least 6 h at hybridization temperature <b>Note:</b> For single copy detection we recommend o/n incubation.

## 2.3 RNA : RNA hybridization

**Hybridization temperature** For RNA:RNA hybridization in general 68° C is the recommended hybridization temperature. The actual hybridization temperature may have to be adjusted depending on the GC content, and homology of probe to target.

**Procedure** In the following table the procedure for a RNA:RNA hybridization is described.  
**Note:** Do not use open trays when working with DIG Easy Hyb buffer.

Step	Action
1	Pre-heat appropriate volume of DIG Easy Hyb (approx. 20 ml/100 cm <sup>2</sup> ) to hybridization temperature.
2	Incubate the membrane for 30 min with gentle agitation. <b>Note:</b> The membrane should be well immersed and covered with DIG Easy Hyb.
3	Denature DIG-labeled RNA probe (100 ng/ml hybridization solution) by boiling for 5 min and rapidly cooling on ice-water.
4	Add to pre-heated DIG Easy Hyb (at least 3.5 ml/100 cm <sup>2</sup> membrane) and mix well but avoid foaming (bubbles may lead to background).
5	Pour off prehybridization solution and immediately add probe/DIG Easy Hyb mixture to membrane. <b>Note:</b> Do not add concentrated probe directly to avoid localized background.
6	Incubate with gentle agitation for at least 6 h at hybridization temperature <b>Note:</b> For detection of rare mRNAs we recommend 16 h incubation time.

## 2.4 DNA : RNA hybridizations

**Hybridization temperature** For DNA:RNA hybridization in general 50°C is the recommended hybridization temperature. The actual hybridization temperature may have to be adjusted depending on the GC content, and homology of probe to target.

**Procedure** In the following table the procedure for a DNA:RNA hybridization is described.  
**Note:** Do not use open trays when working with DIG Easy Hyb buffer.

Step	Action
1	Pre-heat appropriate volume of DIG Easy Hyb (approx. 20 ml/100 cm <sup>2</sup> ) to hybridization temperature.
2	Incubate the blot for 30 min with gentle agitation. <b>Note:</b> The membrane should be well immersed and covered with DIG Easy Hyb.
3	Denature DIG-labeled DNA probe (5-25 ng/ml for DNA-probes, 100 ng/ml for RNA-probes) by boiling for 10 min and rapidly cooling on ice-water.
4	Add to pre-heated DIG Easy Hyb (at least 3.5 ml/100 cm <sup>2</sup> membrane) and mix well but avoid foaming (bubbles may lead to background).
5	Pour off prehybridization solution and immediately add probe/DIG Easy Hyb mixture to membrane. <b>Note:</b> Do not add concentrated probe directly to avoid localized background.
6	Incubate with gentle agitation for at least 6 h at hybridization temperature. <b>Note:</b> For detection of rare mRNAs we recommend 16 h incubation time

## 2.5 Hybridization with DIG labeled oligonucleotide probes

**Hybridization temperature** The hybridization temperature is calculated as follows:  
Calculate  $T_m$  of the oligonucleotide probe by summing up 4° C for each G and C and 2° C for each T or A.  
Perform prehybridization and hybridization at 10° C below evaluated  $T_m$ .

**Multiple locus fingerprinting probes** For multiple locus fingerprinting probes we recommend 2 to 4 h hybridization time. Unspecific competitor DNA like DNA, MB grade from fish sperm (Cat. No. 1 467 140) should be added at a concentration of 50 µg/ml.

**Procedure** In the following table the procedure for hybridization with DIG labeled oligonucleotide probes is described.

**Note:** For tailed oligonucleotides add 0.1 mg/ml poly (A) and 5 µg/ml poly d(A) to the prehybridization and hybridization to prevent unspecific hybridization signals caused by the tails.

Do not use open trays when working with DIG Easy Hyb buffer.

Step	Action
1	Pre-heat appropriate volume of DIG Easy Hyb (approx. 20 ml/ 100 cm <sup>2</sup> ) to hybridization temperature.
2	Incubate the blot for 30 min with gentle agitation. <b>Note:</b> The membrane should be well immersed and covered with DIG Easy Hyb.
3	Hybridize with 0.1-2 pmol tailed oligonucleotide/ml of hybridization solution or 1-10 pmol of end-labeled oligonucleotide. Use at least 3.5 ml DIG Easy Hyb per 100 cm <sup>2</sup> of membrane.
4	Add to pre-heated DIG Easy Hyb (at least 3.5 ml/100 cm <sup>2</sup> membrane) and mix well but avoid foaming (bubbles may lead to background).
5	Pour off prehybridization solution and immediately add probe/DIG Easy Hyb mixture to membrane.
6	Incubate with gentle agitation for 1-6 h at hybridization temperature. <b>Note:</b> For detection of rare mRNAs we recommend 16 h incubation time.

## 2.6 Plaque/ colony hybridization

**Additional reagents required**

- 10% SDS (w/v)
- 20× SSC: 3 M NaCl, 0.3 M sodium citrate, pH 7.0

**Hybridization temperature** The appropriate hybridization temperature is calculated according to G/C content and percent homology of probe to target DNA with the following equation:

$$T_m = 49.82 + 0.41 (\% G+C) - (600/l)$$

[l = length of hybrid in bp]

$$T_{opt.} = T_m - (20^\circ \text{ to } 25^\circ \text{C})$$

The actual hybridization temperature  $T_{opt.}$  with DIG Easy Hyb is 20–25°C below  $T_m$ .

**Procedure**

- The following volumes are calculated for the use of a 275 ml volume roller-bottle.
- The hybridization temperature is given for a 100% homologous probe with 50% G/C content.
- Please make certain that the membranes do not stick to each other and are sufficiently covered with hybridization solution.

**Note:** Do not use open trays when working with DIG Easy Hyb buffer.

Step	Action
1	<ul style="list-style-type: none"> <li>Place 3 membrane discs (82 mm Ø) in a roller bottle and add 60 ml DIG Easy Hyb.</li> <li>Prehybridize for 1 h at 42°C in a hybridization oven for roller bottles.</li> </ul>
2	Denature the labeled probe (25 ng/ml hybridization solution) by boiling for 5 min at 95-100°C and rapidly place on ice.
3	Mix the denatured probe with DIG Easy Hyb, prewarmed to hybridization temperature (5-25 ng/ml).
4	Remove the prehybridization solution and add 6 ml of the probe/DIG Easy Hyb mixture.
5	Incubate for 2 h at 42°C. <b>Note:</b> The hybridization solution with the DIG-labeled probe is stable at -15 to -25°C for more than 12 months and can be reused several times when freshly denatured.

## 2.7 Immunological Detection

The procedure for the immunological detection of DIG-labeled nucleic acids is described in the pack inserts of the DIG Luminescent Detection Kit (Cat. No. 1 363 514) or the DIG Wash and Block Buffer Set (Cat. No. 1 585 762). These pack inserts are available from our website <http://www.roche-applied-science.com>.

## 2.8 Post Hybridization Washes, Stripping and Rehybridization

**Post hybridization washes** Please find in the following table the procedure for the post hybridization washes.

Step	Action
1	Wash 2× 5 min in ample 2× SSC; SDS 0.1% at 15 - 25°C.
2	Wash 2× 15 min in 0.1 × SSC; SDS 0.1% at 68° C under constant agitation.

### Stripping and rehybridization

Please refer to the following table.  
**Note:** When stripping and rehybridization of blots is planned, the membrane should not dry off at any time.

### Caution

Work in a fume hood

Step	Action
1 (only color detection)	Pre-heat dimethylformamide in a waterbath to 50-60°C and incubate the membrane until the color (NBT/BCIP) is washed off. <b>Note:</b> DMF is volatile and can be ignited above 67°C.
2	Rinse membrane briefly in sterile double distilled water.
3	Wash for 2× 20 min in 0.2 M NaOH, SDS, 0.1% (w/v) at 37°C under constant agitation.
4	Equilibrate briefly in 2× SSC.
5	Prehybridize and incubate with second probe.

## 3.2 Reference

1 Itakura, K. et al. (1984) Annu. Rev. Biochem. **53**, 323

\* available from Roche Applied Science

## 3.3 Related products

### Kits

The use of DIG Easy Hyb is recommended in combination with the following DIG kits and replaces the there-mentioned hybridization solutions.

Product	Pack Size	Cat. No
DIG DNA Labeling Kit	40 labeling reactions	1 175 033
DIG RNA Labeling Kit	2 × 10 reactions	1 277 073
DIG Oligonucleotide 3'-End Labeling Kit	25 reactions	1 362 372
DIG Oligonucleotide Tailing Kit	25 reactions	1 417 231
DIG DNA Labeling and Detection Kit	25 labeling reactions and 50 blots (10× 10 cm <sup>2</sup> )	1 093 657
DIG Nucleic Acid Detection Kit	40 blots (10× 10 cm <sup>2</sup> )	1 175 041
DIG Luminescent Detection Kit for Nucleic Acids	50 blots (10× 10 cm <sup>2</sup> )	1 363 514
DIG Wash and Block Buffer Set	1 set	1 585 762

### Single reagents

Product	Pack Size	Cat. No.
Nylon membranes, positively charged	10 sheets (20 × 30 cm)	1 209 272
	20 sheets (10 × 15 cm)	1 209 299
	1 roll (0.3 × 3 m)	1 417 240
Nylon Membranes for Colony and Plaque Hybridization	50 filters (Ø 82 mm)	1 699 075
	50 filters (Ø 132 mm)	1 699 083
DNA, MB-grade	500 mg (50 ml)	1 467 140
Hybridization bags	50 bags	1 666 649

### For your further information:

Roche Applied Science offers a large selection of products for the non-radioactive labeling and detection of nucleic acids.

For a complete overview, please visit and bookmark our "DIG Reagents and Kits for Non-Radioactive Nucleic Acid Labeling and Detection" Special Interest Site at <http://www.roche-applied-science.com/DIG>

## How to contact Roche Applied Science

[www.roche-applied-science.com](http://www.roche-applied-science.com)

to order, solve technical queries, find product information, or contact your local sales representative.

[www.roche-applied-science.com/pack-insert/1603558a.pdf](http://www.roche-applied-science.com/pack-insert/1603558a.pdf)

Please visit our new Online Technical Support Site



Roche Diagnostics GmbH  
 Roche Applied Science  
 Nonnenwald 2  
 82372 Penzberg  
 Germany