

Uni-ZAP[®] XR Library

Custom and Premade Libraries

INSTRUCTION MANUAL

Revision #078001d

BN #937111-12

STORAGE CONDITIONS

Custom Library: 4°C

Premade Library: -80°C

Bacterial Glycerol Stocks: -80°C

Helper Phage: -80°C

For in Vitro Use Only

NOTICE TO PURCHASER

The Uni-ZAP® XR vector is covered by Stratagene's United States Patent No. 5,128,256. The purchase of this vector includes a limited, nonexclusive license under such patent rights to use the vector for the cloning, expression and characterization of genes. This license does not grant rights to (1) use the Uni-ZAP XR vector for the reproduction, amplification or modification of the vector; (2) offer the Uni-ZAP XR vector or any derivative thereof for resale; (3) distribute or transfer the Uni-ZAP XR vector or any derivative thereof to any third party; or (4) incorporate the Uni-ZAP XR vector or any derivative thereof in any genomic or cDNA library for resale, distribution or transfer to any third party. No other license, express, implied or by estoppel, is granted. For information concerning the availability of licenses to reproduce and/or modify the Uni-ZAP XR vector, please contact Stratagene's Technical Services Department at 1-800-424-5444.

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Stratagene. Stratagene shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

United States and Canada

Stratagene Cloning Systems
11011 North Torrey Pines Road
La Jolla, CA 92037
Telephone (619) 535-5400
Order Toll Free (800) 424-5444
Technical Services (800) 424-5444
Internet techservices@stratagene.com
World Wide Web www.stratagene.com

Stratagene European Contacts

Location	Telephone	Fax	Technical Services
Austria	660 312 526	660 312 527	017 956 7036
Belgium	0800 96078	0800 96024	027 13 12 11
Germany	0130 840 911	0130 762 088	0699 509 6197
Netherlands	0800 023 0446	0800 023 0447	0800 023 0448
Switzerland	0800 830 250	0800 825 225	01 800 9045
United Kingdom	0800 585 370	0800 783 0889	0171 365 1056

Distributors

All other countries, please contact your local distributor (see *Stratagene Distributors* in this instruction manual for a complete listing).

Uni-ZAP® XR Library

CONTENTS

Materials Provided.....	1
Additional Materials Required	1
Introduction.....	2
Host Strains	4
Helper Phage	5
Preparing the Host Strains	7
Host Strain Genotypes	7
Growing and Maintaining the Host Strains	7
Preparing a -80°C Bacterial Glycerol Stock	8
Color Selection by IPTG and X-gal.....	8
Preparing Plating Cultures.....	8
Day 1	8
Day 2.....	8
Titering Procedure	9
Amplifying the Library.....	10
Day 1	10
Day 2.....	10
Day 3.....	11
Performing Plaque Lifts.....	11
Titering.....	11
Plating	11
Lifting	12
Hybridizing and Screening	12
Antibody Screening Protocol.....	13
In Vivo Excision of the pBluescript Phagemid from the Uni-ZAP XR Vector	13
In Vivo Excision Protocols Using the ExAssist Helper Phage	13
with SOLR Strain	13
Single-Clone Excision Protocol.....	14
Mass Excision Protocol	15
Troubleshooting.....	16
Preparation of Media And Reagents.....	17
Related Stratagene Products	18
Stratagene Distributors	20
References.....	21
Endnotes	21

Uni-ZAP[®] XR Library

MATERIALS PROVIDED

Materials provided	Quantity	
	Premade Library	Custom Library
Amplified premade library constructed in the Uni-ZAP [®] XR vector ^a	1 ml	—
Custom library constructed in the Uni-ZAP [®] XR vector ^b	—	—
Host strains ^c		
XL1-Blue MRF ⁻ strain	0.5-ml bacterial glycerol stock	0.5-ml bacterial glycerol stock
SOLR [™] strain	0.5-ml bacterial glycerol stock	0.5-ml bacterial glycerol stock
f1 helper phage ^d		
ExAssist [™] interference-resistant helper phage	1 ml	1 ml
VCSM13 Interference-Resistant Helper Phage	1 ml	1 ml

^a Premade libraries have been amplified one time and frozen in the presence of 7% DMSO. Upon arrival, store at -80°C. Do not pass through more than two freeze-thaw cycles.

^b Amplified custom libraries are shipped at 4°C without DMSO. Upon arrival, add DMSO to 7%, aliquot, and store at -80°C. Do not pass through more than two freeze-thaw cycles.

^c Use the SOLR strain for plating excised phagemids and the XL1-Blue MRF⁻ strain for all other manipulations. For host strain storage conditions, see *Preparing the Host Strains*.

^d Retiter after 1 month. (Take care not to contaminate the Uni-ZAP[®] XR vector with this high-titer filamentous helper phage.) Store at -20°C. Stratagene recommends VCSM13 interference-resistant helper phage for single stranded rescue. ExAssist[™] interference-resistant helper phage has α -complementing β -galactosidase sequences which may interfere with sequencing or site-directed mutagenesis where oligonucleotide primers hybridize to β -galactosidase sequences (e.g., M13-20 primer).

ADDITIONAL MATERIALS REQUIRED

Isopropyl-1-thio- β -D-galactopyranoside (IPTG)
5-Bromo-4-chloro-3-indoyl- β -D-galactopyranoside (X-gal)
Nitrocellulose membrane

Revision #078001d

Copyright © 1998 by Stratagene.

INTRODUCTION

The Uni-ZAP[®] XR vector system combines the high efficiency of lambda library construction and the convenience of a plasmid system with blue–white color selection. The Uni-ZAP XR vector (see Figure 1) is double digested with *EcoR* I and *Xho* I and will accommodate DNA inserts from 0 to 10 kb in length. The Uni-ZAP XR vector can be screened with either DNA probes or antibody probes and allows in vivo excision of the pBluescript[®] phagemid (see Figure 2), allowing the insert to be characterized in a plasmid system¹. The polylinker of the pBluescript phagemid has 21 unique cloning sites flanked by T3 and T7 promoters and a choice of 6 different primer sites for DNA sequencing. The phagemid has the bacteriophage f1 origin of replication, allowing rescue of single-stranded DNA, which can be used for DNA sequencing or site-directed mutagenesis. Unidirectional deletions can be made with exonuclease III and mung bean nuclease by taking advantage of the unique positioning of 5' and 3' restriction sites. Transcripts made from the T3 and T7 promoters generate riboprobes useful in Southern and Northern blotting, and the *lacZ* promoter may be used to drive expression of fusion proteins suitable for Western blot analysis or protein purification.






FIGURE 1. Map of the Uni-ZAP[®] XR insertion vector.

The library was synthesized using the ZAP-cDNA[®] synthesis method.² The linker–primer was designed with a GAGA sequence to protect the *Xho* I restriction enzyme recognition site and an 18-base poly(dT) sequence. The restriction site allows the finished cDNA to be inserted into the vector unidirectionally in the sense orientation with respect to the *lacZ* promoter. The linker–primer is a 50-base oligonucleotide with the following sequence:



The adaptors are comprised of 9- and 13-mer oligonucleotides, which are complementary to each other and have an *Eco*R I cohesive end. The adaptors have the following sequence:



The pBluescript SK(–) plasmid in the Uni-ZAP XR vector contains the N-terminus of the *lacZ* gene, which can be α -complemented by the specific host strain used. There are 36 amino acids from the MET sequence to the *Eco*R I site. A total of 131 amino acids are coded for, but this is interrupted by the large polylinker.

Host Strains

The RecA[–] *E. coli* host strain XL1-Blue MRF' is supplied with the Uni-Zap XR library.³ Because the Uni-ZAP XR vector does not require a *supF* genotype, the amplified library grows very efficiently on the XL1-Blue MRF' strain. In addition, use of the correct host strain is important when working with the Uni-ZAP XR vector as the F' episome present in the XL1-Blue MRF' strain serves three purposes.

First, the $\Delta M15$ *lacZ* gene present on the F' episome is required for the β -galactosidase-based nonrecombinant selection strategy. When cDNA is present in the polylinker, expression from the *lacZ* gene is disrupted and white plaques are produced. In contrast, without insert in the polylinker, the amino terminus of β -galactosidase is expressed and nonrecombinants can be scored visually by the presence of blue plaques. To produce an enzymatically active β -galactosidase protein, two domains are required: the α -region expressed by the vector and the $\Delta M15$ *lacZ* domain expressed by the F' episome. These two domains fold to form a functional protein, the α -region complementing the missing amino acids resulting from the $\Delta M15$ mutation. Therefore, in order to utilize the nonrecombinant selection strategy, the correct host strain must be used to produce a functional β -galactosidase protein.

Second, the F' episome expresses the genes forming the F' pili found on the surface of the bacteria. Without pili formation, filamentous phage (i.e., M13 or f1) infection could not occur. Because the conversion of a recombinant Uni-ZAP XR clone to a pBluescript phagemid requires superinfection with a filamentous helper phage, the F' episome is required for in vivo excision (see *In Vivo Excision of the pBluescript Phagemid from the Uni-ZAP XR Vector*).

Third, the F' episome contains the *lac* repressor (*lacI^q* gene), which blocks transcription from the *lacZ* promoter in the absence of the inducer isopropyl-1-thio- β -D-galactopyranoside (IPTG) (see *Related Stratagene Products*). This repressor is important for controlling expression of fusion proteins which may be toxic to the *E. coli*. Because the presence of the *lacI^q* repressor in the *E. coli* host strain can potentially increase the representation or completeness of the library, XL1-Blue MRF' is useful for screening the amplified library.

Note *The strains used for the Lambda gt11 vector (i.e., Y1088, Y1089, and Y1090) are not suitable for use with the Uni-ZAP XR vector because these strains contain the plasmid pMC9, a pBR322 derivative, which contains many of the same sequences as those found in the phagemid portion of the Uni-ZAP XR vector. Using these strains with the Uni-ZAP XR vector could result in recombination between the homologous sequences.*

Helper Phage

Two different helper phages are provided with the Uni-Zap XR library: (1) the ExAssist interference-resistant helper phage with SOLR strain⁴ and (2) the VCSM13 helper phage. The ExAssist interference-resistant helper phage with SOLR strain is designed to allow efficient in vivo excision of the pBluescript phagemid from the Uni-ZAP XR vector while preventing the problems that can be associated with helper phage co-infection. The ExAssist helper phage contains an amber mutation that prevents replication of the phage genome in a nonsuppressing *E. coli* strain (e.g., SOLR cells). Only the excised phagemid can replicate in the host, removing the possibility of co-infection from the ExAssist helper phage. The ExAssist helper phage cannot be used for single-stranded rescue due to its inability to replicate in the SOLR strain. The other helper phage, VCSM13 helper phage, is recommended for single-stranded rescue procedures from the excised pBluescript phagemids.

Storing the Helper Phage

The ExAssist helper phage and the VCSM13 helper phage may be stored at -80° , -20° , or 4° C. Some loss of titer may occur after 2 months of storage at 4° C. The stability will increase to 1 year at -80° C. If the titer drops over time, prepare a fresh high-titer stock of the helper phage as outlined in *Amplification*.

Titering the Helper Phage

If the ExAssist helper phage has been stored at 4° C for >1 month or passed through a freeze-thaw cycle, titer the helper phage with XL1-Blue MRF' cells.

1. Dilute the phage (10^{-4} – 10^{-7}) in TE buffer and combine 1 μ l of each dilution with 200 μ l of XL1-Blue MRF' cells ($OD_{600} = 1.0$) as outlined in steps 2–5 below. Expect titers of $<10^{10}$ pfu/ml.
2. Incubate the helper phage and the XL1-Blue MRF' cells for 15 minutes at 37° C to allow the phage to attach to the cells.
3. Add 3 ml of NZY top agar,[§] melted and cooled to $\sim 48^{\circ}$ C, and then pour immediately onto prewarmed NZY agar plates.[§]

4. Incubate the plates overnight at 37°C.

Note *ExAssist* plaques will have a cloudier appearance than lambda phage plaques.

5. To determine the titer [in plaque-forming units per milliliter (pfu/ml)], use the following formula:

$$\left[\frac{\text{Number of plaques (pfu)} \times \text{dilution factor}}{\text{Volume plated } (\mu\text{l})} \right] \times 1000 \mu\text{l} / \text{ml}$$

where the volume plated (in microliters) refers to the volume of the helper phage solution added to the cells.

Amplifying the Helper Phage

1. Transfer a colony of XL1-Blue or XL1-Blue MRF' cells from a fresh LB-tetracycline plate into 10 ml of 2× YT broth (see *Preparation of Media and Reagents*) in a 50-ml conical tube.
2. Incubate the conical tube with shaking at 37°C until growth reaches an OD₆₀₀ of 0.3.

Note *An OD₆₀₀ of 0.3 corresponds to 2.5 × 10⁸ cells/ml.*

3. Add the helper phage at a multiplicity of infection (MOI) of 20:1 (phage-to-cells ratio).

Note *When amplifying VCSM13 helper phage, add kanamycin to a final concentration of 25 μg/ml after 30 minutes of growth.*

4. Grow the culture at 37°C for 8 hours.
5. Heat the conical tube at 65°C for 15 minutes.
6. Spin down the cell debris and transfer the supernatant to a fresh conical tube.
7. The titer of the supernatant should be between 7.5 × 10¹⁰ and 1.0 × 10¹² pfu/ml for ExAssist helper phage or between 1.0 × 10¹¹ and 1.0 × 10¹² pfu/ml for VCSM13 helper phage.
8. Add dimethylsulfoxide (DMSO) to a final concentration of 7% (v/v) and store at –80°C.

For further details about helper phage titering or amplification, please see *Titering Procedure* or reference 5.

PREPARING THE HOST STRAINS

Host Strain Genotypes

Host strains	Genotype
SOLR strain ^a	e14 ⁻ (McrA ⁻) Δ(mcrCB-hsdSMR-mrr)171 sbcC recB recJ uvrC umuC::Tn5 (Kan ^r) lac gyrA96 relA1 thi-1 endA1 λ ^R [F' proAB lac ^r ZΔM15] Su ⁻ (nonsuppressing)
XL1-Blue MRF' strain	Δ(mcrA)183 Δ(mcrCB-hsdSMR-mrr)173 endA1 supE44 thi-1 recA1 gyrA96 relA1 lac [F' proAB lac ^r ZΔM15 Tn10 (Tet ^r)]

^a Use the SOLR strain for excision only.

Growing and Maintaining the Host Strains

The bacterial host strains are shipped as bacterial glycerol stocks. For the appropriate media, please refer to the following table:

Note *The mcrA, mcrCB and mrr mutations prevent restriction of methylated DNA, making XL1-Blue MRF' compatible with cloning both genomic DNA and cDNA constructed using Stratagene's ZAP-cDNA[®] synthesis kit.*

Host strain	Agar plates for bacterial streak	Medium for bacterial glycerol stock	Medium for bacterial cultures for titering phage (final concentration)
SOLR strain ^a	LB-kanamycin ^b	LB-kanamycin ^b	LB with 0.2% (w/v) maltose-10 mM MgSO ₄
XL1-Blue MRF' strain	LB-tetracycline ^b	LB-tetracycline ^b	LB with 0.2% (w/v) maltose-10 mM MgSO ₄

^a Use the SOLR strain for excision only.

^b See *Preparation of Media and Reagents*.

On arrival, prepare the following from the bacterial glycerol stock using the appropriate media as indicated in the previous table:

Note *The host strains may thaw during shipment. The vials should be stored immediately at -20° or -80°C, but most strains remain viable longer if stored at -80°C. It is also best to avoid repeated thawing of the host strains in order to maintain extended viability.*

1. Revive the stored cells by scraping off splinters of solid ice with a sterile wire loop.
2. Streak the splinters onto an LB agar plate (see *Preparation of Media and Reagents*) containing the appropriate antibiotic, if one is necessary.
3. Incubate the plate overnight at 37°C.
4. Seal the plate with Parafilm[®] laboratory film and store the plate at 4°C for up to 1 week.
5. Restreak the cells onto a fresh plate every week.

Preparing a -80°C Bacterial Glycerol Stock

1. In a sterile 50-ml conical tube, inoculate 10 ml of appropriate liquid medium containing antibiotic with one colony from the plate. Grow the cells to late log phase.
2. Add 4.5 ml of a sterile glycerol–liquid medium solution (5 ml of glycerol + 5 ml of appropriate medium) to the bacterial culture from step 1. Mix well.
3. Aliquot into sterile centrifuge tubes (1 ml/tube).

This preparation may be stored at -20°C for 1–2 years or at -80°C for more than 2 years.

Color Selection by IPTG and X-gal

The color selection by α -complementation with the Uni-ZAP XR vector requires higher amounts of IPTG[‡] and X-gal[‡] for generation of the blue color. Transcription and translation of the fusion protein are normal, but the large polylinker present within the pBluescript phagemid, which is present in the Uni-ZAP XR vector, is partly responsible for the reduced activity of the β -galactosidase protein—not the promoter. As would be expected, the copy number of the Uni-ZAP XR vector is much less per cell than the copy number of pBluescript phagemids. However, it is important to note that the color assay is used only for determining the ratio of recombinants to nonrecombinants within a newly constructed library and is not used for any other manipulations.

PREPARING PLATING CULTURES

Day 1

1. Inoculate 50 ml of LB broth supplemented with 0.2% (w/v) maltose and 10 mM MgSO_4 in a sterile flask with a single colony of the appropriate bacterial host.

Note *Do not add antibiotic to the overnight culture or to the titering plates.*

2. Grow overnight with shaking at 30°C . This temperature ensures that the cells will not overgrow. Phage can adhere to dead cells as well as to live ones and can lower the titer.

Day 2

3. Spin the cells down in a sterile conical tube for 10 minutes at 2000 rpm.
4. Carefully decant the media off the cell pellet and **gently** resuspend the pellet in ~15 ml of 10 mM MgSO_4 . (Do not vortex.)
5. Dilute the cells to $\text{OD}_{600} = 0.5$ with 10 mM MgSO_4 . Approximately 600 μl of $\text{OD}_{600} = 0.5$ cells are needed for each 150-mm plate and 200 μl of $\text{OD}_{600} = 0.5$ cells for each 100-mm plate.
6. Use freshly prepared cells for titering and amplifying the library.

[‡] See *Related Stratagene Products*

TITERING PROCEDURE

A background test can be completed by plating several hundred plaques on a plate [see *Color Selection by IPTG and X-gal*]. Add 15 μl of 0.5 M IPTG (in water) and 50 μl of 250 mg/ml X-gal [in dimethylformamide (DMF)] to 2–3 ml of NZY top agar, melted and cooled to $\sim 48^\circ\text{C}$. The higher concentrations of IPTG and X-gal used in the plating often result in the formation of a precipitate, which disappears after incubation. Stratagene suggests that the IPTG and X-gal should be added separately, with mixing in between additions, to the NZY agar plates. Plaques are visible after incubation for 6–8 hours at 37°C . Nonrecombinant (background) plaques are blue, while recombinant plaques are white.

1. Prepare the host bacteria as outlined in *Preparing Plating Cultures*.
2. Make dilutions of the lambda phage in SM buffer.
3. Add 1 μl of the lambda phage to 200 μl of host cells diluted in 10 mM MgSO_4 to $\text{OD}_{600} = 0.5$. If desired, also add 1 μl of a 1:10 dilution of the packaged material in SM buffer to 200 μl of host cells.
4. Incubate the phage and bacteria for 15 minutes at 37°C to allow the phage to attach to the cells.
5. Add the following components:
 - 2–3 ml of NZY top agar (melted and cooled to $\sim 48^\circ\text{C}$).
 - 15 μl of 0.5M IPTG (in water)
 - 50 μl of X-gal [250 mg/ml (in DMF)]
6. Plate immediately onto NZY agar plates and allow the plates to set for 10 minutes. Invert the plates and incubate at 37°C .
7. Plaques should be visible after 6–8 hours, although color detection requires overnight incubation. Background plaques are blue and should be $< 1 \times 10^5$ pfu/ μg of arms, while recombinant plaques will be white (clear) and should be 10–100-fold above the background.

Note *Primary libraries can be unstable; therefore, amplification of the libraries is recommended immediately.*

AMPLIFYING THE LIBRARY

It is usually desirable to amplify libraries prepared in lambda vectors to make a large, stable quantity of a high-titer stock of the library. However, more than one round of amplification is not recommended, since slower growing clones may be significantly underrepresented.

Note *The premade library has been through one round of amplification.*

The following protocol is recommended for amplifying the Uni-ZAP XR library:

Day 1

1. Prepare the host strains as outlined in *Preparing Plating Cultures*

Note *In order to obtain the highest amplification efficiency, use freshly prepared host strains.*

Day 2

2. Dilute the cells to an OD₆₀₀ of 0.5 in 10 mM MgSO₄. Use 600 µl of cells at an OD₆₀₀ of 0.5/150-mm plate.
3. Combine aliquots of the packaged mixture or library suspension containing $\sim 5 \times 10^4$ pfu of bacteriophage with 600 µl of host cells at an OD₆₀₀ of 0.5 in Falcon® 2059 polypropylene tubes. To amplify 1×10^6 plaques, use a total of 20 aliquots (each aliquot contains 5×10^4 plaques/150-mm plate).

Note *Do not add more than 300 µl of phage/600 µl of cells.*

4. Incubate the tubes containing the phage and host cells for 15 minutes at 37°C.
5. Mix 6.5 ml of NZY top agar, melted and cooled to $\sim 48^\circ\text{C}$, with each aliquot of infected bacteria and spread evenly onto a freshly poured 150-mm NZY agar plate.
6. Incubate the plates at 37°C for 6–8 hours. Do not allow the plaques to get larger than 1–2 mm. On completion, the plaques should be touching.
7. Overlay the plates with ~ 8 –10 ml of SM buffer. Store the plates at 4°C overnight (with *gentle* rocking if possible). This allows the phage to diffuse into the SM buffer.

Day 3

8. Recover the bacteriophage suspension from each plate and pool it into a sterile polypropylene container. Rinse the plates with an additional 2 ml of SM buffer and pool. Add chloroform to a 5% (v/v) final concentration. Mix well and incubate for 15 minutes at room temperature.
9. Remove the cell debris by centrifugation for 10 minutes at $500 \times g$.
10. Recover the supernatant and transfer it to a sterile polypropylene container. If the supernatant appears cloudy or has a high amount of cell debris, repeat steps 8 and 9. If the supernatant is clear, add chloroform to a 0.3% (v/v) final concentration and store at 4°C . Stratagene recommends storing aliquots of the amplified library in 7% (v/v) DMSO at -80°C .
11. Check the titer of the amplified library using host cells and serial dilutions of the library. (Assume $\sim 10^9$ – 10^{11} pfu/ml.)

PERFORMING PLAQUE LIFTS

Titering

Titer the library using fresh host cells to determine the concentration as outlined in *Plating and Titering*.

Plating

Plate the library at 50,000 pfu/plate on a series of large 150-mm NZY agar plates, which are at least 2 days old, following the procedure outlined below.

1. Add the equivalent of 50,000 pfu/plate to 600 μl of host cells at an OD_{600} of 0.5.
2. Incubate the bacteria and phage mixture at 37°C for 15 minutes to allow the phage to attach to the cells.
3. Add 6.5 ml of NZY top agar to the bacteria and phage mixture.
4. Quickly pour the plating culture onto an NZY agar plate and carefully swirl the plate to distribute the cells evenly. (Use 20 plates to screen 1×10^6 pfu.)
5. Invert the plates and incubate at 37°C for ~ 8 hours.
6. Chill the plates for 2 hours at 4°C to prevent the NZY top agar from sticking to the nitrocellulose membrane.

Lifting

Note Use forceps and wear gloves for the following steps.

1. Place a nitrocellulose membrane[‡] onto each NZY agar plate for 2 minutes to allow the transfer of the phage particles to the membrane. Use a needle to prick through the membrane and agar for orientation. (If desired, waterproof ink in a syringe needle may be used.)

Notes If making duplicate nitrocellulose membranes, allow the second membrane to transfer for ~4 minutes.

Pyrex® dishes are convenient for the following steps. All solutions should be at room temperature.

- a. Denature the nitrocellulose membrane after lifting by submerging the membrane in a 1.5 M NaCl and 0.5 M NaOH denaturation solution for 2 minutes.

Note If using charged nylon, wash with gloved fingertips to remove the excess NZY top agar.

- b. Neutralize the nitrocellulose membrane for 5 minutes by submerging the membrane in a 1.5 M NaCl and 0.5 M Tris-HCl[‡] (pH 8.0) neutralization solution.
 - c. Rinse the nitrocellulose membrane for no more than 30 seconds by submerging the membrane in a 0.2 M Tris-HCl (pH 7.5) and 2× SSC buffer (see *Preparation of Media and Reagents*) solution.
2. Blot briefly on a Whatman® 3MM paper.
 3. Crosslink the DNA to the membranes using the autocrosslink setting on the Stratalinker® UV crosslinker[‡] (120,000 μJ of UV energy) for ~30 seconds. Alternatively, oven bake at 80°C for ~1.5–2 hours.
 4. Store the agar plates at 4°C.

HYBRIDIZING AND SCREENING

Following the preparation of the membranes for hybridization, perform prehybridization, probe preparation, hybridization, and washes for either oligonucleotide probes or double-stranded probes and then expose the membranes to film as outlined in standard methodology texts.^{5,6} Following these procedures, perform secondary and tertiary screenings also as outlined in the standard methodology texts.^{5,6} After an isolate is obtained, refer to Sambrook *et al.*⁵ for suggested phage miniprep and majorprep procedures.

[‡] See *Related Stratagene Products*.

ANTIBODY SCREENING PROTOCOL

A complete manual for immunoscreening is supplied with Stratagene's *picoBlue*[™] immunoscreening kit. This kit is available with goat anti-rabbit antibodies and goat anti-mouse antibodies (see *Related Stratagene Products*).

In Vivo EXCISION OF THE pBLUESCRIPT PHAGEMID FROM THE UNI-ZAP XR VECTOR

The Uni-ZAP XR vector is designed to allow simple, efficient in vivo excision and recircularization of any cloned insert contained within the lambda vector to form a phagemid containing the cloned insert. This in vivo excision depends on the placement of the DNA sequences within the lambda phage genome and on the presence of a variety of proteins, including f1 bacteriophage-derived proteins. The f1 phage proteins recognize a region of DNA normally serving as the f1 bacteriophage origin of replication. This origin of replication can be divided into two overlying parts: (1) the site of initiation and (2) the site of termination for DNA synthesis.⁷ These two regions are subcloned separately into the Uni-ZAP XR vector. The lambda phage (target) is made accessible to the f1-derived proteins by simultaneously infecting a strain of *E. coli* with both the lambda vector and the f1 bacteriophage.

Inside *E. coli*, the "helper" proteins (i.e., proteins from f1 or M13 phage) recognize the initiator DNA that is within the lambda vector. One of these proteins then nicks one of the two DNA strands. At the site of this nick, new DNA synthesis begins and duplicates whatever DNA exists in the lambda vector "downstream" (3') of the nicking site. DNA synthesis of a new single strand of DNA continues through the cloned insert until a termination signal, positioned 3' of the initiator signal, is encountered within the constructed lambda vector. The single-stranded DNA molecule is circularized by the gene II product from the f1 phage, forming a circular DNA molecule containing the DNA between the initiator and terminator. In the case of the Uni-ZAP XR vector, this includes all sequences of the pBluescript SK(-) phagemid and the insert, if one is present. This conversion is the "subcloning" step, since all sequences associated with normal lambda vectors are positioned outside of the initiator and terminator signals and are not contained within the circularized DNA. In addition, the circularizing of the DNA automatically recreates a functional f1 origin as found in f1 bacteriophage or phagemids.

In Vivo EXCISION PROTOCOLS USING THE EXASSIST HELPER PHAGE WITH SOLR STRAIN

The ExAssist helper phage with SOLR strain is designed to allow efficient excision of the pBluescript phagemid from the Uni-ZAP XR vector, while eliminating problems associated with helper phage co-infection. The ExAssist helper phage contains an amber mutation that prevents replication of the phage genome in a nonsuppressing *E. coli* strain such as SOLR cells. This allows only the excised phagemid to replicate in the host, removing the possibility of co-infection from the ExAssist helper phage. Since the ExAssist helper phage cannot replicate in the SOLR strain, single-stranded rescue cannot be performed in this strain using this helper phage.

Mass excision can be used to generate subtraction libraries and subtracted DNA probes.

Single-Clone Excision Protocol

Day 1

1. Core the plaque of interest from the agar plate and transfer the plaque to a sterile microcentrifuge tube containing 500 μ l of SM buffer and 20 μ l of chloroform. Vortex the microcentrifuge tube to release the phage particles into the SM buffer. Incubate the microcentrifuge tube for 1–2 hours at room temperature or overnight at 4°C. (This phage stock is stable for up to 6 months at 4°C.)
2. Grow separate overnight cultures of XL1-Blue MRF' and SOLR cells in LB broth (see *Preparation of Media and Reagents*), supplemented with 0.2% (w/v) maltose and 10 mM MgSO₄, at 30°C.

Day 2

3. Gently spin down the XL1-Blue MRF' and SOLR cells (1000 \times g). Resuspend the XL1-Blue MRF' and SOLR cells at an OD₆₀₀ of 1.0 in 10 mM MgSO₄.
4. Combine the following components in a Falcon 2059 polypropylene tube:

200 μ l of XL1-Blue MRF' cells at an OD₆₀₀ of 1.0
250 μ l of phage stock (containing $>1 \times 10^5$ phage particles)
1 μ l of the ExAssist helper phage ($>1 \times 10^6$ pfu/ μ l)

5. Incubate the Falcon 2059 polypropylene tube at 37°C for 15 minutes.
6. Add 3 ml of LB broth and incubate the Falcon 2059 polypropylene tube for 2.5–3 hours at 37°C with shaking. Because clonal representation is not relevant, single-clone excision reactions can be safely performed overnight.

Note *The turbidity of the media is not indicative of the success of the excision.*

7. Heat the Falcon 2059 polypropylene tube at 65–70°C for 20 minutes and then spin the tube at 1000 \times g for 15 minutes.
8. Decant the supernatant into a sterile Falcon 2059 polypropylene tube. This stock contains the excised pBluescript phagemid packaged as filamentous phage particles. (This stock may be stored at 4°C for 1–2 months.)
9. To plate the excised phagemids, add 200 μ l of freshly grown SOLR cells from step 3 (OD₆₀₀ = 1.0) to two 1.5-ml microcentrifuge tubes. Add 100 μ l of the phage supernatant from step 8 above to one microcentrifuge tube and 10 μ l of the phage supernatant to the other microcentrifuge tube.
10. Incubate the microcentrifuge tubes at 37°C for 15 minutes.
11. Plate 200 μ l of the cell mixture from each microcentrifuge tube on LB–ampicillin agar plates (50 μ g/ml) (see *Preparation of Media and Reagents*) and incubate the plates overnight at 37°C.

Due to the high-efficiency of the excision process, it may be necessary to titrate the supernatant to achieve single-colony isolation.

Colonies appearing on the plate contain the pBluescript double-stranded phagemid with the cloned DNA insert. Helper phage will not grow, since helper phage is unable to replicate in the Su⁻ (nonsuppressing) SOLR strain and does not contain ampicillin-resistance genes. SOLR cells are also resistant to lambda phage infection, thus preventing lambda phage contamination after excision.

To maintain the pBluescript phagemid, streak the colony on a new LB–ampicillin agar plate. For long-term storage, prepare a bacterial glycerol stock and store at -80°C.

VCSM13 helper phage is recommended for the the single-stranded rescue procedure. The single-stranded rescue procedure can be found in Stratagene's *pBluescript[®] Exo/Mung DNA Sequencing System Instruction Manual*.

Mass Excision Protocol

Day 1

1. Grow separate overnight cultures of XL1-Blue MRF['] and SOLR cells in LB broth, supplemented with 0.2% (w/v) maltose and 10 mM MgSO₄, at 30°C.

Day 2

2. Gently spin down the XL1-Blue-MRF['] and SOLR cells (1000 × g). Resuspend the XL1-Blue MRF['] and the SOLR cells in 10 mM MgSO₄ to an OD₆₀₀ of 1.0 (8 × 10⁸ cells/ml).
3. In a 50-ml conical tube, combine a portion of the amplified lambda bacteriophage library with XL1-Blue MRF['] cells at a MOI of 1:10 lambda phage-to-cell ratio. Excise 10- to 100-fold more lambda phage than the size of the primary library to ensure statistical representation of the excised clones. Add ExAssist helper phage at a 10:1 helper phage-to-cells ratio to ensure that every cell is co-infected with lambda phage and helper phage.

For example, use

10⁷ pfu of the lambda phage (i.e., 10- to 100-fold above the primary library size)
10⁸ XL1-Blue MRF['] cells (1:10 lambda phage-to-cell ratio, noting that an OD₆₀₀
of 0.3 corresponds to 2.5 × 10⁸ cells/ml)
10⁹ pfu of ExAssist helper phage (10:1 helper phage-to-cells ratio)

4. Incubate the conical tube at 37°C for 15 minutes to allow absorption.
5. Add 20 ml of LB broth and incubate the conical tube for 2.5–3 hours at 37°C with shaking.

Notes *Incubation times for mass excision in excess of 3 hours may alter the clonal representation.*

The turbidity of the media is not indicative of the success of the excision.

6. Heat the conical tube at 65–70°C for 20 minutes.
7. Spin down the debris at 1000 × g for 10 minutes and then decant the supernatant into a sterile conical tube.
8. To titer the excised phagemids, combine 1 µl of this supernatant with 200 µl of SOLR cells from step 2 in a 1.5-ml microcentrifuge tube.
9. Incubate the microcentrifuge tube at 37°C for 15 minutes.
10. Plate 100 µl of the cell mixture onto LB–ampicillin agar plates (100 µg/ml) and incubate the plates overnight at 37°C.

Note *It may be necessary to further dilute the cell mixture to achieve single-colony isolation.*

At this stage, colonies may be selected for plasmid preps, or the cell mixture may be plated directly onto filters for colony screening.

TROUBLESHOOTING

Observations	Possible causes	Suggestions
The number of colonies is too low	The molar ratios of lambda phage to cells to helper phage is critical	Verify that the titer on the tubes is current and correct and use only calibrated pipettor
	Excision efficiencies are directly related to the Uni-ZAP XR phage titer	If an excision is unsuccessful, prepare a high-titer stock of the phage and repeat the excision procedure
	Poor rescue may be a result of toxic cDNA clones which can be isolated in lambda vectors but not in plasmid vectors	The ABLE C strain ^a and the ABLE™ K strain ^a reduce the copy number of common cloning vectors by ~4- and 20-fold, respectively, enhancing the probability that a toxic clone will be propagated. Positive clones observed on initial screening as lambda plaques can be excised and introduced into the ABLE strains. Excised phagemid libraries can also be screened directly in the ABLE strains
	Chloroform, which is added after packaging to prevent bacterial contamination, lyses the <i>E. coli</i> before the helper phage can infect and excise	Be sure to spin down the chloroform completely prior to removing an aliquot of the lambda phage for in vivo excision

^a See *Related Stratagene Products*

PREPARATION OF MEDIA AND REAGENTS

Note All media must be autoclaved before use.

LB Agar (per Liter)

10 g of NaCl
10 g of tryptone
5 g of yeast extract
20 g of agar
Add deionized H₂O to a final volume of 1 liter
Adjust pH to 7.0 with 5 N NaOH
Autoclave
Pour into petri dishes (~25 ml/100-mm plate)

LB Broth

10 g of NaCl
10 g of tryptone
5 g of yeast extract
Add deionized H₂O to a final volume of 1 liter
Adjust pH to 7.0 with 5 N NaOH
Autoclave

LB-Tetracycline Broth (per Liter)

Prepare 1 liter of LB broth
Autoclave
Cool to 55°C
Add 1.5 ml of 10 mg/ml tetracycline (filter-sterilized)
Store broth in a dark, cool place as tetracycline is light-sensitive

NZY Broth (per Liter)

5 g of NaCl
2 g of MgSO₄ · 7H₂O
5 g of yeast extract
10 g of NZ amine (casein hydrolysate)
Add deionized H₂O to a final volume of 1 liter
Adjust the pH to 7.5 with NaOH
Autoclave

LB-Ampicillin Agar (per Liter)

Prepare 1 liter of LB agar
Autoclave
Cool to 55°C
Add 10 ml of 10 mg/ml ampicillin (filter-sterilized)
Pour into petri dishes (~25 ml/100-mm plate)

LB-Kanamycin Broth (per Liter)

Prepare 1 liter of LB broth
Autoclave
Cool to 55°C
Add 7.5 ml of 10 mg/ml kanamycin (filter-sterilized)

LB-Tetracycline Agar (per Liter)

Prepare 1 liter of LB agar
Autoclave
Cool to 55°C
Add 1.5 ml of 10 mg/ml tetracycline (filter-sterilized)
Pour into petri dishes (~25 ml/100-mm plate)
Store plates in a dark, cool place or cover plates with foil if left out at room temperature for extended time periods as tetracycline is light-sensitive

NZY Agar (per Liter)

5 g of NaCl
2 g of MgSO₄ · 7H₂O
5 g of yeast extract
10 g of NZ amine (casein hydrolysate)
15 g of agar
Add deionized H₂O to a final volume of 1 liter
Adjust the pH to 7.5 with NaOH
Autoclave
Pour into petri dishes (~80 ml/150-mm plate)

LB Top Agar (per Liter)

Prepare 1 liter of LB broth
Add 0.7% (w/v) agarose
Autoclave

NZY Top Agar (per Liter)

Prepare 1 liter of NZY broth
Add 0.7% (w/v) agarose
Autoclave

SM Buffer (per Liter)

5.8 g of NaCl
2.0 g of MgSO₄ · 7H₂O
50.0 ml of 1 M Tris-HCl (pH 7.5)
5.0 ml of 2% (w/v) gelatin
Add deionized H₂O to a final volume of 1 liter
Autoclave

20× SSC Buffer (per Liter)

175.3 g of NaCl
88.2 g of sodium citrate
800.0 ml of deionized H₂O
Adjust to pH 7.0 with a few drops of 10 N NaOH
Add deionized H₂O to a final volume of 1 liter

2× YT Broth (per Liter)

10 g of NaCl
10 g of yeast extract
16 g of tryptone
Add deionized H₂O to a final volume of 1 liter
Adjust to pH 7.5 with NaOH
Autoclave

RELATED STRATAGENE PRODUCTS**Escherichia coli Strains**

ABLE™ C strain and ABLE™ K strain [Catalog #200305 (2 × 0.5-ml glycerol stock)]
ABLE™ C strain [Catalog #200306 (0.5-ml glycerol stock)]
ABLE™ K strain [Catalog #200307 (0.5-ml glycerol stock)]
SOLR™ strain [Catalog #200298 (0.5-ml glycerol stock)]
XL1-Blue MRF' strain [Catalog #200301 (0.5-ml glycerol stock)]

Electroporation-Competent Cells

Epicurian Coli® XL1-Blue MRF' electroporation-competent cells (Catalog #200158)
Epicurian Coli® ABLE™ electroporation-competent cell kit (Catalog #200160)
Epicurian Coli® ABLE™ C electroporation-competent cells (Catalog #200161)
Epicurian Coli® ABLE™ K electroporation-competent cells (Catalog #200162)

Competent Cells

Epicurian Coli® ABLE™ competent cell kit (Catalog #200170)
Epicurian Coli® ABLE™ C competent cells (Catalog #200171)
Epicurian Coli® ABLE™ K competent cells (Catalog #200172)
Epicurian Coli® XL1-Blue MRF' supercompetent cells (Catalog #200230)

Equipment and Labware

Stratalinker® UV crosslinker [Catalog #400071 (1800) and #400075 (2400)]
Duralon-UV™ membranes (Catalog #420100–#420105)
Nitrocellulose Membranes (Catalog #420106–#420108)

Systems and Kits

Prime-It® II random primer labeling kit (Catalog #300385)

KinAce-It™ kinasing kit (Catalog #200390)

picoBlue™ antibody screening kit [Catalog #200371 (goat anti-rabbit) and #200372 (goat anti-mouse)]

QuikHyb® hybridization solution [Catalog #201220 (250 ml) and #201221 (1 liter)]

DNA Ligation Kit (Catalog #203003)

Fine Chemicals

IPTG [Catalog #300127 (1 g)]

X-gal [Catalog #300200 (250 mg) and #300201 (1 g)]

Salmon Sperm DNA (Catalog #201190)

STRATAGENE DISTRIBUTORS

Country	Distributor name	Telephone number	Fax number
Argentina	Tecnolab S.A.	541 555 0010	541 553 3331
Australia	Integrated Sciences Toll free	(02) 9417 7866 (800) 252204	(02) 9417 5066
Austria	See inside front cover		
Brazil	Instrucom	11 5561 1771	11 530 0895
Canada	See inside front cover		
Denmark	AH Diagnostics	86 10 10 55	86 16 15 33
Egypt	Clinilab	202 3518763	202 3781507
Finland	Kebo Finland	358 9 804 551	358 9 804 55200
France	Ozyme	1 34 60 24 24	1 34 60 92 12
Germany	See inside front cover		
Hong Kong	Line Analytics Ltd.	2578 5839	2807 2674
India	Wipro Ltd.	91 11 3325677	91 11 3738675
Ireland	B.M. Brownes Ltd.	353 1 295 3401	353 1 295 3818
Israel	Getter Photo Supplies (2000) Ltd.	3 576 1555	3 752 3620
Italy	Eppendorf s.r.l.	02-58.01.34.09	02-58.01.34.38
Japan	Funakoshi Co., Ltd. Toyobo Co., Ltd. (Tokyo) Toyobo Co., Ltd. (Osaka)—Main office	3 5684 1622 3 3660 4819 6 348 3785/6 348 3788	3 5684 1633 3 3660 4887 6 348 3322
Korea	Koram Biotech Corp.	2 556 0311	2 556 0828
Mexico	Bioselec	355 8928	556 6943
New Zealand	LabSupply Pierce (NZ) Ltd.	9 443 5867	9 444 7314
Norway	MedProbe A.S.	47 22 20 01 37	47 22 20 01 89
Portugal	Biocontec	1 361 3620	1 362 5615
Republic of China	Merck Taiwan Ltd.	886 2 2521 9331	886 2 2536 7734
Singapore	ITS Science & Medical Pte. Ltd.	273 0898	273 0810
South Africa	Whitehead Scientific Supplies	21 981 1560	21 981 5789
Spain	Cultek	91 729 03 33	91 358 17 61
Sweden	AH Diagnostics AB Toll free	86 80 08 45 800 10299	86 80 04 35
Switzerland	See inside front cover		
Thailand	ITS (Thailand) Co., Ltd.	2 308 0611	2 308 0612
United Kingdom	See inside front cover		
United States	See inside front cover		

REFERENCES

1. Short, J. M., and Sorge, J. A. (1992) *Meth. Enzymol.* 216:495–508.
2. Short, J. M., *et al.* (1988) *Nucleic Acids Res.* 16:7583–7600.
3. Bullock, W. O., Fernandez, J. M., and Short, J. M. (1987) *Biotechniques* 5:376–379.
4. Hay, B., and Short, J. M. (1992) *Strategies* 5(1): 16–18.
5. Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989) "Molecular Cloning: A Laboratory Manual," Second Ed. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York.
6. Ausubel, F.M., Brent, R., Kingston, R.E. Moore, D.D., Seidman, J.G., Smith, J.A., and Struhl, K. (1990) "Current Protocols in Molecular Biology." Wiley, New York.
7. Dotto, G. P., Horiuchi, K., and Zinder, N. D. (1984) *J. Mol. Biol.* 172: 507–521.

ENDNOTES

Uni-ZAP[®], pBluescript[®], Epicurian Coli[®], Stratalinker[®], Prime-It[®], QuikHyb[®], and ZAP-cDNA[®] are registered trademarks of Stratagene in the United States.

SOLR, Duralon-UV, ExAssist, KinAce-It, picoBlue, and ABLE are trademarks of Stratagene.

Falcon[®] is a registered trademark of Becton Dickinson Labware.

GenBank[®] is a registered trademark of the National Institutes of Health.

Parafilm[®] is a registered trademark of American Can Company.

Pyrex[®] is a registered trademark of Corning Glass Works.

Whatman[®] is a trademark of Whatman Paper Ltd.