

Please note the following changes to the BD In-Fusion Dry-Down Kit User Manual (PT3754-1): Add components to Cat. No. 639605 that are shown in bold below.

II. List of Components

Cat. No. 639605

- 96 rxn BD In-Fusion™ Dry-Down Mix
- **12** **Optically Clear PCR Cap Strips**
- **1** **Microseal A Film**
- 5 µl pDNR-Dual (100 ng/µl, linear)
- 5 µl pDNR-CMV (100 ng/µl, linear)
- 5 µl 1.1 kb Control Insert (25 ng/µl)